

SafeCube Dispenser™

**Reduce The Risk Of Accidental Drips, Drops & Spills!
& Enhances Regulatory Compliance!**

**The Safest method for filling smaller containers
from five gallon cubitainers of laboratory chemicals
eg., 10% Formalin, Disinfectants, & Reagents.**

System Benefits

- No more Dripping Spigots all over the floor from a 5-gallon cube hanging over the edge of the counter.
- No need to balance a container while filling, simply place container to be filled in filling tray and turn on spigot.
- Safe while filling containers and then simply close it up for safe storage when your not filling.
- Control the Drips, Drops & Accidental Spills with Secondary Spill Containment. (Main Cabinet holds five (5) Gallons)
- Liquid Detection Alarm to immediately signal working personnel of any accidental spills.



Perfect For Use In:

- Surgery and Labor & Delivery
- Ambulatory Surgery
- Histo-Pathology Laboratory
- Clinical and Research Laboratory
- College & University Laboratory
- Morgue & Autopsy Suites
- Housekeeping & Central Services

PRODUCT # SCD-100A

Measurements: 20" deep x 18 1/2" wide x 17" high