

**PDI - PROFESSIONAL DISPOSABLES, INC  
HEALTHCARE DIVISION OF  
NICE-PAK PRODUCTS, INC**

**TWO NICE-PAK PARK  
ORANGEBURG NY 10962-1376  
845-365-1700**

REVISION DATE: 10/03

**I. PRODUCT IDENTIFICATION:** SANI-DEX TOWELETTES - D35400, Q35372

PRODUCT/TRADE NAME:

HAZARD RATING(NFPA)

HEALTH: 1

FLAMMABILITY: 3

REACTIVITY: 0

SPECIFIC: NONE

EMERGENCY OR INFORMATION TELEPHONE NO: 845-365-1700 (M-F DAYTIME) AT OTHER TIMES, CONTACT THE LOCAL POISON CONTROL CENTER

CHEMICAL NAME: PROPRIETARY MIXTURE

**II HAZARDOUS INGREDIENTS PER 29 CFR 1910.1200**

Hazardous Ingredients	%	ACGIH TLV	CAS NUMBER
SD ALCOHOL 40	40	1000 PPM	64-17-5
PARACHLOROMETAXYLENOL 0.5		NOT ESTABLISHED	88-04-0

**III-PHYSICAL/CHEMICAL CHARACTERISTICS**

**III-PHYSICAL/CHEMICAL CHARACTERISTICS**

COLOR/ODOR/APPEARANCE: TOWELETTE SATURATED WITH

CLEAR LIQUID WITH FLORAL/ALCOHOL ODOR..

BOILING POINT: N/A

FLASH POINT: 80<sub>o</sub>F

VAPOR DENSITY N/A

EVAPORATION RATE: N/A

SOLUBILITY IN WATER: COMPLETE

SPECIFIC GRAVITY (H<sub>2</sub>O = 1): 0.940

**IV FIRE & EXPLOSION HAZARD DATA**

FLASH POINT(Method Used): 80<sup>o</sup>F TAG CLOSED CUP : LEL: 2.0 UEL: 12.0

EXTINGUISHING MEDIA: DRY CHEMICAL OR ALCOHOL TYPE FOAM, CARBON DIOXIDE

SPECIAL FIRE FIGHTING PROCEDURES: HANDLE AS FLAMMABLE LIQUID. USE RESPIRATORY PROTECTION FOR FIRE FIGHTING PERSONNEL

UNUSUAL FIRE AND EXPLOSION HAZARDS: CLASS 3 FLAMMABILITY

**V-REACTIVITY DATA**

STABILITY: STABLE

CONDITIONS TO AVOID NONE

INCOMPATIBILITY: NONE

HAZARDOUS DECOMPOSITION OR BYPRODUCT: NONE

POLYMERIZATION: WILL NOT OCCUR

CONDITIONS TO AVOID: NONE

SANI-DEX

## **VI-HEALTH HAZARD DATA**

EFFECTS OF OVER EXPOSURE:

SKIN IF RASH APPEARS, DISCONTINUE USE

EYES : WILL CAUSE EYE STING IF SPLASHED

INHALATION: NONE

INGESTION NOT A NORMAL ROUTE OF ENTRY

EMERGENCY AND FIRST AID PROCEDURES:

SKIN CONTACT: DISCONTINUE USE IF RASH OR IRRITATION OCCURS

EYE CONTACT: FLUSH WITH COLD WATER FOR 15 MINUTES

INHALATION: REMOVE TO FRESH AIR

INGESTION: IF SWALLOWED, CALL POISON CONTROL CENTER

TARGET ORGANS; N/A

## **VII-SPILL AND DISPOSAL PROCEDURE**

SPILL CONTROL: ELIMINATE ALL SOURCES OF IGNITION

WASTE DISPOSAL METHOD: FLUSH SPILLS WITH WATER. HANDLE AS FLAMMABLE LIQUID. FOLLOW LOCAL, STATE AND FEDERAL REGULATIONS.

HANDLING AND STORAGE: STORE AWAY FROM HEAT AND SOURCES OF IGNITION

## **VIII-CONTROL MEASURES/PROTECTION**

RESPIRATION: USE IF PERMISSABLE EXPOSURE LEVEL IS EXCEEDED WHEN HANDLING BULK LIQUID

VENTILATION: RECOMMENDED

PROTECTIVE GLOVES: NONE REQUIRED

EYE PROTECTION: ONLY IF SPLASHING IS EXPECTED

HYGIENIC PRACTICES:GOOD HOUSEKEEPING PRACTICES SHOULD BE FOLLOWED.

OTHER: DON'T ALLOW LARGE QUANTITIES OF WASTE TO ACCUMULATE

## **IX TRANSPORT/SHIPPING**

DOT SHIPPING NAME: CONSUMER COMMODITY

TECHNICAL SHIPPING NAME: N/A

DOT SHIPPING CLASSIFICATION: "ORM-D"

DOT ID NO.: N/A

DOT LABEL REQUIREMENTS: N/A

UN/NA NUMBER REGULATIONS: N/A

REPORTABLE QUANTITY: N/A

## **X DISCLAIMER**

THE INFORMATION FURNISHED HEREIN IS BELIEVED TO BE ACCURATE AND REPRESENTS THE BEST DATA CURRENTLY AVAILABLE TO US. NO WARRANTY, EXPRESSED OR IMPLIED IS MADE AND NICE-PAK PRODUCTS, INC. ASSUMES NO LEGAL RESPONSIBILITY OR LIABILITY RESULTING FROM ITS USE.